

## Components




CoolSpeK main body	1	<b>Standard Accessories</b>
Temperature Controller	1	1 Silicon Tube ( ID 5mm x OD 9mm x 3 Meters )
Liquid Nitrogen Reservoir	1	Set of Tubes for Gas Flow (urethane tube, connector of taper pipe threads, flow valve)
*Cuvette is not supplied		Tool Kit
		User's Manual
<b>Outer Size of the main body (without options)</b>		
146.5mm(H) x 90mm(W) x 111mm(D)		

## Specifications

Liquid Nitrogen Reservoir	Stainless, 2L
Low-Temperature Sample Chamber	Aluminium, polyurethane foam for thermal insulation
Optical Windows	Quartz, 3-way
Suitable Cell	Outer dimension 12.5mm x 12.5mm
Temperature Control	Regulated liq. N <sub>2</sub> flow
Temperature Range	-80°C to +100°C
Volume of Liq. N <sub>2</sub> Consumption	1L / hour
Precision of Temperature Control	±1°C or ±0.5% of indication value, whichever is greater (error of the sensor not included)
Quantity of Dew Condensation	Less than 0.10D / hour at -80°C with UNISOKU spectrophotometer
Temperature Sensor	Resistance thermometer sensor (Pt-100 Class B)
Functions	2 built-in heaters in the main body One for prevention of dew condensation on optical windows. Another for temperature control.
Cryogen Used	Liquid Nitrogen


## Application Examples

### Time-dependent UV-Vis Spectra



Spectral changes for the reaction of O<sub>2</sub> and a Cu (I) complex at -80°C. Formation of a m-peroxo dinuclear copper (II) complex can be easily monitored by the time-dependent UV-Vis spectrum at low temperature. (Courtesy of Prof. Shinobu Itoh of Osaka University)

### Variable-temperature UV-Vis, CD and Fluorescence Spectra



Valuable temperature UV-Vis, circular dichroism, and fluorescence spectra of (R)-1-methylheptyl mellitate (0.03mM) in pentane at 25°C to -110°C. FL: Excitation at 270nm.

(Courtesy of Prof. Yoshihisa Inoue of Osaka University and Prof. Takehiko Wada of Tohoku University)

\*Instrument components are subject to change without prior notice for improvement in performance.

# Cryostat for Spectrophotometer

# CoolSpeK UV/CD

## USP-203 Series



For UV-Vis spectrophotometers  
CoolSpeK UV USP-203

For Circular Dichroism (CD) Spectrophotometers  
CoolSpeK CD USP-203CD

User friendly

-80°C to +100°C in units of 0.1°C

Striking performance with liquid solution

Kinetics

Light weight and compact footprint

CoolSpeK allows you to obtain spectra and chemical kinetic rates under low-temperature condition easily.

# Features

CoolSpeK is attached to a sample compartment of commercial UV-Vis or fluorescence spectrophotometers by using an adaptor. Various kinds of options are prepared to meet customers requests.

## Chemical Kinetics

CoolSpeK is able to measure a chemical reaction process by injecting a reagent and using a magnetic stirrer (option).

## User friendly

CoolSpeK does not require vacuum pumps, and can be used in atmospheric pressure. It cools the sample to low temperature by flowing liquid N<sub>2</sub> from a reservoir, and controls the temperature precisely by auto-regulating the flow.

## Low Dew Condensation

Special structure and heating function prevent dew condensation: Less than 0.05OD / hour at -80°C with UNISOKU spectrophotometer, after heat treatment.

## Temperature Range

-80°C ~ Room Temperature ~ +100°C

## Extensive Customer Base

Shipped over 500 sets globally (2022)

## Light Weight and Compact Footprint


CoolSpeK is compact and can be placed inside the sample compartments of commercial UV-Vis or fluorescence spectrophotometers by using adapters. Cuvettes with cap for spectrophotometers can be used for this cryostat.

\*The lid of sample compartments in some spectrophotometers may not close when the cryostat is attached.

## Variety of Options

There are plenty of options. Magnetic stirrer, various base adaptors, cuvette adaptor, solid sample holder, and so on.

## Example of Use



## Installation Example



## OPTIONS

### Magnetic Stirrer (No. CS-AT-SM)

This is attached at the bottom of CoolSpeK main body. The spin speed is controlled by the controller.  
\*This magnetic stirrer cannot be attached to some spectrophotometers.

### Standard Cuvette (No. CS-CL-U1)

Four side transparent quartz cuvette with a screw cap (light path length: 10mm)

### Magnetic stirring Bar (No. CS-SB-F27)

This is used in the Standard Cuvette (No. CS-CL-U1). Specified Tolerance of Teflon® Walls  
High Strength Magnetization (φ 2mm x 7mm)

### Adaptor for 1mm Light Path Length Cell (No. CS-CL-H1) Adaptor for 2mm Light Path Length Cell (No. CS-CL-H2)

This adaptor fills a gap between the cuvette and the cryostat internal body. A spring in the adaptor makes thermal conduction better.  
\*Does not include cuvette

### Adaptor for 1mm Light Path Length Cuvette for Fluorescence (No. CS-CL-F1)

This adaptor enables you to acquire the surface fluorescence of samples with high absorption in 1mm light path length cuvette.  
\*Does not include cuvette

### Solid Sample Holder for Transmittance (No. CS-KT-H00-32 or H07-42)

Suitable sample: solid sheet or powder  
Suitable dimension: φ 10mm or 10mm x 10mm, Thickness: 0 ~ 3.2mm or 0.7 ~ 4.2mm  
\*Does not include cuvette  
\*Please contact us if your sample dimension differs from above.

### Solid Sample Holder for Fluorescence (No. CS-KF-H00-32 or H07-42)

Suitable sample: solid sheet or powder  
Suitable dimension: φ 13mm x 13mm, Thickness: 0 ~ 3.2mm or 0.7 ~ 4.2mm  
\*Please contact us if your sample dimension differs from above.

## Systems Interfaceable with CoolSpeK UV

UNISOKU RSP-1000/TSP-1000  
Agilent Technologies Agilent (Cary) 8453 / 8454  
Agilent Technologies Cary 50/60  
Agilent Technologies Cary 3500 Flexible Module  
(Cary 3500 Compact / Multicell not available)  
Agilent Technologies Cary 4000 / 5000 / 6000  
Beckman DU-7400  
JASCO V-500/600/700 Series  
JASCO FP-6000/8000 Series  
JASCO J-720/820/1500  
Hitachi U-2000/3000/4000 Series  
Hitachi F-4500/7000

Horiba FluoroMax/Log  
Horiba FluoroCube  
Perkin Elmer Lambda Series  
Shimadzu UV1800/2000/3000 series  
Shimadzu RF-5300/6000  
Sinco S-3100  
Edinburgh Instruments FLS920 / FLS1000  
Edinburgh Instruments LP980(TAS)  
※ FLS adaptor "product code F-5376" manufactured by EI separately required

\*Some spectrophotometers might not be able to shut out the room light completely when the cryostat is attached.  
\*Please contact us if your spectrophotometer is not listed here.