# 2021 UNISOKU NEWSLETTER

Let the global  
pandemic  
be blown out!

# Business Policy

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UNISOKU Co., Ltd. was founded by Dr. Toshihiko Nagamura in 1974 with the spirit of his frontier mentality, "Provide unique measurement systems to the society". In 2010, we joined Tokyo Instruments(TII) Group, the worldwide trading group for advanced scientific instruments. Under the former President, Shoji Suruga, UNISOKU promoted globalization and capital investment and succeeded in expanding our sales worldwide.

We keep listening to our customers, continuously improve our systems, and visit them until their systems work properly and positive outcomes are obtained. This is our basic policy to build a bond with our customers. We strive to develop unique instruments that can provide significant impacts and create new values in the society, and continue to grow in the constantly changing world. A variety of innovative measurement systems have been continuously produced from UNISOKU, including not only very low temperature/high magnetic field STM systems but also tip-enhanced Raman spectroscopy systems possessing molecule-scale spatial resolution, cryogen-free STM systems, and optical instruments such as picoTAS, which enables us to measure over a wide range of time scales including "Gap Time Region" by utilizing our own RIPT method.

It is our mission and the origin of our growth to listen to our customers and provide them happiness and invaluable experience by offering solutions to their problems with our skills and knowledge. We continue to challenge and never give up until our products meet the customers' needs, and strive to further improve our products and service.

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President and CEO Yutaka Miyatake

## Our Policy for Delivery and Repair Visit Work due to COVID-19

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Due to the pandemic COVID-19, we sincerely apologize for any late deliveries and suspending our on-site visit work. Regarding restarting our on-site visit work, UNISOKU must follow the immigration policy announced by the Government of Japan for each country. Although it is difficult to restart traveling abroad depending on area, we do our best to keep this inconvenience to a minimum by offering online support and/or local support by the SPECS-TII group. We deeply appreciate your kind understanding in this regard. We are also grateful to our suppliers for their efforts to provide us products and keep their production lines running during this difficult situation.

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## SPECS-TII Group

### SPECS-TIIグループ

UNISOKU joined TII group as a subsidiary company in 2010. In 2017, TII Group and SPECS GmbH established a joint company, SPECS-TII GmbH in Switzerland, for strengthening their worldwide sales and support. SPECS-TII GmbH then established local firms in China, Russia and the United States for sales and support in those areas. UNISOKU and SPECS-TII aim to provide further reassurance and service for our customers using our global network.

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## Corporate Philosophy

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UNISOKU contributes to the development of science and technology by providing customers with measurement systems that meet their exploring minds.

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## The Code of Conduct

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#### Challenge courageously

Even in difficult situations, we try something new with our strong will.

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#### Think positively and look from others' standpoints

We express our gratitude. We achieve goals through mutual cooperation. We think from the standpoints of customers, partners and colleagues, and provide value for them.

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#### Create values

To meet ever-changing social needs, we aim to create values that provide happiness and impressive experience to customers by our knowledge and skills.

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

#### Pursue happiness

We work to share happiness with our customers and colleagues.

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2月  
FEB

2nd fluorescence lifetime option for picoTAS was delivered to Osaka University.  
picoTAS用蛍光寿命オプション2台納品 46期2/2台目(阪大)

3月  
MAR

The RIPT method development in JST-SENTAN program received an S grade (the highest evaluation).  
JST先端計測(RIPT法開発)でS評価獲得

UNISOKU delivered Room Temperature THz-SNOM system to Fudan University (through TII).  
復旦大学(TII経由)常温THz-SNOMシステムを出荷

4月  
APR

UNISOKU started teleworking for the first time due to COVID-19.  
コロナウィルス感染防止対策で初めてテレワークを導入

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

## Production Division

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## Design Division

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We decided to hire him as a full-time employee from a contractual worker because we recognize his sincere attitude and speed to learn a new job. We expect him to have more experience in production and contribute to the production of high quality systems.

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His experience in mechanical design is attractive for us. We expect him to understand mechanical structures of UNISOKU products and become a designer who can propose appropriate designs to our customers in the future.

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
## Product Development News in 2020

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## Highly Hydrogen-Sensitive Thermal Desorption Spectroscopy (TDS) System

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~10<sup>3</sup> higher sensitivity to hydrogen than other commercial TDS systems

Arbitrary heating rate between 10-100°C/min for 150-1000°C



This work was in collaboration with Prof. Hosono's group (Tokyo Institute of Technology) and supported by Ichimura Foundation for New Technology.

## Cryogen-Free STM (USM1800)

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
STM image of HOPG

- Base temperature below  $T_{\text{sample}} = 5 \text{ K}$
- Low vibration noise level realized by our original design

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- Kelvin probe by KP Technology Ltd.
- Sample transfer & exchange in UHV conditions
- Variable sample temperature

## Applications


In-situ work function & surface potential measurements for thin films and clean surfaces prepared in vacuum

Distributor : Tokyo Instruments, Inc.

## Raman Spectroscopy System for Thin Films Under Various Environmental Conditions

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Where to deliver : Japan Atomic Energy Agency Dr. Terasawa Tomo-o

Combination with Nanofinder® by Tokyo Instruments, Inc.



- Sample transfer & exchange in UHV conditions
- Raman spectroscopy measurements under gaseous/heating environments

## Applications

- Thin films grown by CVD/MBE
- Variable temperature measurements
- Surface-adsorbed molecules at low temperatures

## Hall Effect Measurement System for the Autonomous Synthesis System

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
Designed for automatic Hall effect measurements and sample transfer by a robot arm




This work was in collaboration with Prof. Hitosugi's group (Tokyo Institute of Technology) and MI-6 Ltd. as a part of their NEDO project.

## Scanning Noise Microscope (SNOiM)

GaAs nano-device

Electron temperature mapping

Real-space imaging of electron temperature with nanometer-scale spatial resolution  
Q. Weng *et al.* Science **360**, 775 (2018).

This work is in collaboration with Dr.s Weng and Komiyama (Kim group in RIKEN).

# Introduction of Products / Users / Publications

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## Publication Stats in 2020

Total number of publications using UNISOKU systems\* = 246 (240 in 2019)

Total impact factors ~ 1898 (1720 in 2019)

Corresponding to 44 Nature papers (40 in 2019)


Impact factor per employee ~ 40 (~36 in 2019)

Approaching to the impact factor of Nature (~43)

The detailed information about the publication list is available on our website.

\*including preprints

Popular Research Fields	Num. of Publications	Average Impact Factor
Molecules	32	7.97
Topological Materials (Majorana)	31	15.77
Graphene (Twisted Bi-layer Graphene)	25	8.48
Transition Metal Dichalcogenides (TMDs)	21	10.69
Low Dimensional Structures (Nanowire, Low Dimensional Superconductivity, etc)	21	7.28
RIPT Transient Absorption	13	8.17
Monatomic Films excluding TMDs	10	10.18
Fe-based Superconductors (FeSe/STO)	9	12.50

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Ultra High Vacuum Low Temperature SPM System

**Solution for long-term measurements with low cost cryogen**

**Application fields:** molecules, transition metal dichalcogenides, low-dimensional structures

## USM1200 Publication List in 2020 (Selected) 論文リスト

**“Reconfigurable Photo-Induced Doping of Two-Dimensional Van Der Waals Semiconductors Using Different Photon Energies”**

S. Seo *et al.*, Nat. Electron. **4**, 38 (2021).

**“Synthesis of Monolayer Blue Phosphorus Enabled by Silicon Intercalation”**

J. L. Zhang *et al.*, ACS Nano **14**, 3687 (2020).

**“Steering Effect of Bromine on Intermolecular Dehydrogenation Coupling of Poly(p-phenylene) on Cu(111)”**

Y. Lin *et al.*, ACS Nano **14**, 17134 (2020).

**“Orbital-Selective Dirac Fermions and Extremely Flat Bands in Frustrated Kagome-Lattice Metal CoSn”**

Z. Liu *et al.*, Nat. Commun. **11**, 4002 (2020).

**“On-surface Synthesis of Nitrogen Substituted Gold-Phosphorus Porous Network”**

J. Zhang *et al.*, Chem. Mater. **32**, 8561 (2020).

**“Molecular-Scale Investigation of the Thermal and Chemical Stability of Monolayer PTCDA on Cu(111) and Cu(110)”**

C. Gu *et al.*, ACS Appl. Mater. Interfaces **12**, 22327 (2020).

**“Single-Molecule Imaging of Dinitrogen Molecule Adsorption on Individual Iron Phthalocyanine”**

C. Gu *et al.*, Nano Res. **13**, 2393 (2020).

**“Realization of a Buckled Antimonene Monolayer on Ag(111) via Surface Engineering”**

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
Ultra High Vacuum Low Temperature SPM System

**Ideal system for combining SPM with optical options and extending to a multi-probe platform**

**Application fields:**

graphene, tip-enhanced Raman spectroscopy, scanning tunneling potentiometry, etc.

## USM1400 Publication List in 2020 (Selected) 論文リスト

**“Sub-Nanometre Resolution in Single-Molecule Photoluminescence Imaging”**

B. Yang *et al.*, Nat. Photonics **14**, 693 (2020).

**“Reversible Superconductor–Insulator Transition in (Li, Fe)OHFeSe Flakes Visualized by Gate-Tunable Scanning Tunneling Spectroscopy”**

R. Yin *et al.*, ACS Nano **14**, 7513 (2020).

**“Experimental Realization of Two-Dimensional Buckled Lieb Lattice”**

H. Feng *et al.*, Nano Lett. **20**, 2537 (2020).

**“Dramatic Enhancement of Tip-Enhanced Raman Scattering Mediated by Atomic Point Contact Formation”**

S. Liu *et al.*, Nano. Lett. **20**, 5879 (2020).

**“Interfacial Superconductivity in FeSe Ultrathin Films on SrTiO<sub>3</sub> Probed by In Situ Independently Driven Four-Point-Probe Measurements”**

A. K. Pedersen *et al.*, Phys. Rev. Lett. **124**, 227002 (2020).

**“Tunable Lattice Reconstruction, Triangular Network of Chiral One-Dimensional States, and Bandwidth of Flat Bands in Magic Angle Twisted Bilayer Graphene”**

Y. Liu *et al.*, Phys. Rev. Lett. **125**, 236102 (2020).

**“The Structure and Mechanism of Large-Scale Indium-Intercalated Graphene Transferred from SiC Buffer Layer”**

T. Hu *et al.*, Carbon **171**, 829 (2021).

**“Phase-Resolved Detection of Ultrabroadband THz Pulses Inside a Scanning Tunneling Microscope Junction”**

M. Muller *et al.*, ACS Photonics **7**, 2046 (2020).

**“Tunable Magnetism of a Single-Carbon Vacancy in Graphene”**

Y. Zhang *et al.*, Sci. Bull. **65**, 194 (2020).


**“Anisotropic Gapping of Topological Weyl Rings in the Charge-Density-Wave Superconductor In<sub>x</sub>TaSe<sub>2</sub>”**

Y. Li *et al.*, Sci. Bull. DOI: 10.1016/j.scib.2020.09.007

## “Sub-Nanometre Resolution in Single-Molecule Photoluminescence Imaging”

B. Yang *et al.*, Nat. Photonics, **14**, 693 (2020). Product used: USM1400

Sub-nanometer resolution imaging with scanning near-field optical microscopy (SNOM) has long been a challenge. Yang *et al.* (Dong group, Univ. of Science and Technology of China) demonstrated spatially and spectrally resolved photoluminescence imaging of a single zinc phthalocyanine (ZnPc) molecule with a spatial resolution down to 0.8 nm. The unprecedented spatial resolution was achieved by (1) a strongly localized-and-enhanced plasmon field using a Ag tip with an atomistic protrusion on the apex, (2) a three-monolayer-thick NaCl dielectric spacer to electronically decouple the ZnPc molecule from a plasmonic substrate Ag(100), and (3) the precise control of the tip-sample tunneling junction with sub-nanometer scale. They revealed intriguing photophysics of molecular fluorescence enhancement and quenching, and subtle plasmon-molecule interactions at the sub-molecular level, providing a new way for optical imaging, spectroscopy and engineering of light-matter interactions at the sub-nanometer level.

(a) Experimental set-up.


(b) Tip-enhanced photoluminescence image of a single ZnPc molecule on NaCl/Ag(100).

(c) Line-profile along AB in (b), showing the spatial resolution of 0.8 nm.

## Best-selling very low temperature SPM system

**Application fields:** topological matters (Majorana), twisted bi-layer graphene, monolayer superconductors etc.


## USM1300 Publication List in 2020 (Selected) 論文リスト

- “Chiral Superconductivity in Heavy-Fermion Metal  $\text{UTe}_2$ ”  
L. Jiao *et al.*, Nature **579**, 523 (2020).
- “Quantum-Limit Chern Topological Magnetism in  $\text{TbMn}_5\text{Sn}_6$ ”  
J. Yin *et al.*, Nature **583**, 533 (2020).
- “Topological Superconductivity in a Van Der Waals Heterostructure”  
S. Kezilebieke *et al.*, Nature **588**, 424 (2020).
- “Evidence for Dispersing 1D Majorana Channels in an Iron-Based superconductor”  
Z. Wang *et al.*, Science **367**, 104 (2020).
- “Nearly Quantized Conductance Plateau of Vortex Zero Mode in an Iron-Based Superconductor”  
S. Zhu *et al.*, Science **367**, 189 (2020).
- “Visualizing Coexisting Surface States in the Weak and Crystalline Topological Insulator  $\text{Bi}_2\text{TeI}$ ”  
N. Avraham *et al.*, Nat. Mater. **19**, 610 (2020).
- “Topological Frustration Induces Unconventional Magnetism in a Nanographene”  
S. Mishra *et al.*, Nat. Nanotechnol. **15**, 22 (2020).
- “Superhigh Uniform Magnetic Cr Substitution in a 2D  $\text{Mo}_2\text{C}$  Superconductor for a Macroscopic-Scale Kondo Effect”  
C. Xu *et al.*, Adv. Mater. DOI: 10.1002/adma.202002825
- “Epitaxial Growth of Single-Phase 1T'- $\text{WSe}_2$  Monolayer with Assistance of Enhanced Interface Interaction”  
W. Chen *et al.*, Adv. Mater. DOI: 10.1002/adma.202004930
- “Interfacial Hydrogen-Bonding Dynamics in Surface-Facilitated Dehydrogenation of Water on  $\text{TiO}_2(110)$ ”  
S. Tan *et al.*, J. Am. Chem. Soc. **142**, 826 (2020).
- “Proximity-Induced Superconductivity in Monolayer  $\text{MoS}_2$ ”  
D. J. Trainer *et al.*, ACS Nano **14**, 2718 (2020).
- “Direct Observation of One-Dimensional Peierls-type Charge Density Wave in Twin Boundaries of Monolayer  $\text{MoTe}_2$ ”  
L. Wang *et al.*, ACS Nano **14**, 8299 (2020).
- “An Anomalous Magneto-Optic Effect in Epitaxial Indium Selenide Layers”  
W. Fu *et al.*, Nano Lett. **20**, 5330 (2020).
- “Relativistic Artificial Molecules Realized by Two Coupled Graphene Quantum Dots”  
Z. Fu *et al.*, Nano Lett. **20**, 6738 (2020).

## “Nearly Quantized Conductance Plateau of Vortex Zero Mode in an Iron-Based Superconductor”

S. Zhu, L. Kong, L. Cao, H. Chen *et al.*, Science **367**, 189 (2020). Product used: USM1300

The search for topological superconductors has been recently one of central issues in the condensed matter physics. In a vortex core of the topological superconductor, Majorana bound states are predicted to exist at zero energy and obey Non-Abelian statistics, offering potential applications in topological quantum computation. However, experimental demonstrations have remained challenging. Zhu *et al.* (Gao group, Chinese Academy of Sciences) performed variable-tunnel-coupled STS measurements at  $T_{\text{eff}} = 377$  mK on vortex cores in  $\text{FeTe}_{0.55}\text{Se}_{0.45}$ , a strong candidate for the topological superconductor. They observed nearly quantized conductance plateaus as a function of tunnel coupling for zero-energy bound states, whereas no plateaus were observed on finite energy bound states originating from conventional quasiparticles. The plateau behavior of the zero-bias conductance is theoretically supported as spectroscopic evidence for the existence of Majorana bound states. Their study moves one step further toward the braiding operation in topological quantum computation.

(a) A variable tunnel coupling STM/STS for a vortex core of  $\text{FeTe}_{0.55}\text{Se}_{0.45}$ .  
(b)  $dI/dV$  spectra taken at the vortex core as a function of the tunnel coupling ( $T_{\text{eff}} = 377$  mK).  
(c) (top) A color-scale plot of (b). (middle) Line-cut at  $V = 0$  meV, showing a conductance plateau.  
(bottom) Line-cuts at high  $V$  values, showing the absence of the conductance plateau.


## USM1300 Publication List in 2020 (Selected) 論文リスト

- “Precise Tuning of Band Structures and Electron Correlations by van der Waals Stacking of One-Dimensional  $\text{W}_6\text{Te}_6$  Wires”  
J. Deng *et al.*, Nano Lett. **20**, 8866 (2020).
- “Mottness Versus Unit-Cell Doubling as the Driver of the Insulating State in 1T-Ta $\text{S}_2$ ”  
C. J. Butler *et al.*, Nat. Commun. **11**, 2477 (2020).
- “Molecular Molds for Regularizing Kondo States at Atom/Metal Interfaces”  
X. Li *et al.*, Nat. Commun. **11**, 2566 (2020).
- “Fermion-Boson Many-Body Interplay in a Frustrated Kagome Paramagnet”  
J. Yin *et al.*, Nat. Commun. **11**, 4003 (2020).
- “Spin-Orbit Quantum Impurity in a Topological Magnet”  
J. Yin *et al.*, Nat. Commun. **11**, 4415 (2020).
- “Long-Range Focusing of Magnetic Bound States in Superconducting Lanthanum”  
H. Kim *et al.*, Nat. Commun. **11**, 4573 (2020).
- “Imaging the Coupling Between Itinerant Electrons and Localised Moments in the Centrosymmetric Skyrmion Magnet  $\text{GdRu}_2\text{Si}_2$ ”  
Y. Yasui *et al.*, Nat. Commun. **11**, 5925 (2020).
- “Single Particle Tunneling Spectrum of Superconducting  $\text{Nd}_{1-x}\text{Sr}_x\text{NiO}_2$  Thin Films”  
Q. Gu *et al.*, Nat. Commun. **11**, 6027 (2020).
- “A New Majorana Platform in an Fe-As Bilayer Superconductor”  
W. Liu *et al.*, Nat. Commun. **11**, 5688 (2020).
- “Atomically Thin Superconductors”  
Z. Li *et al.*, Small DOI: 10.1002/smll.201904788
- “Valley Polarization and Inversion in Strained Graphene via Pseudo-Landau Levels, Valley Splitting of Real Landau Levels, and Confined States”  
S. Li *et al.*, Phys. Rev. Lett. **124**, 106802 (2020).
- “Movable Valley Switch Driven by Berry Phase in Bilayer-Graphene Resonators”  
Y. Liu *et al.*, Phys. Rev. Lett. **124**, 166801 (2020).
- “Many-Body Resonance in a Correlated Topological Kagome Antiferromagnet”  
S. Zhang *et al.*, Phys. Rev. Lett. **125**, 046401 (2020).
- “Strain Tunable Semimetal–Topological-Insulator Transition in Monolayer 1T'- $\text{WTe}_2$ ”  
C. Zhao *et al.*, Phys. Rev. Lett. **125**, 046801 (2020).
- “Direct Visualization of Ambipolar Mott Transition in Cuprate  $\text{CuO}_2$  Planes”  
Y. Zhong *et al.*, Phys. Rev. Lett. **125**, 077002 (2020).
- “Local Berry Phase Signatures of Bilayer Graphene in Intervalley Quantum Interference”  
Y. Zhang *et al.*, Phys. Rev. Lett. **125**, 116804 (2020).
- “Multiple In-Gap States Induced by Topological Surface States in the Superconducting Topological Crystalline Insulator Heterostructure  $\text{Sn}_{1-x}\text{Pb}_x\text{Te-Pb}$ ”  
H. Yang *et al.*, Phys. Rev. Lett. **125**, 136802 (2020).
- “Superconductivity in a Hole-Doped Mott-Insulating Triangular Adatom Layer on a Silicon Surface”  
X. Wu *et al.*, Phys. Rev. Lett. **125**, 117001 (2020).
- “Observation of Discrete Conventional Caroli–de Gennes–Matricon States in the Vortex Core of Single-Layer  $\text{FeSe/SrTiO}_3$ ”  
C. Chen *et al.*, Phys. Rev. Lett. **124**, 097001 (2020).

“Chiral Superconductivity in Heavy Fermion Metal  $\text{UTe}_2$ ”

L. Jiao *et al.*, Nature **579**, 523 (2020). Product used: USM1300

Spin triplet superconductors consisting of electron pairs with spin 1 and an odd-parity wavefunction are topologically nontrivial and expected to provide an ideal platform for realizing Majorana edge states but extremely rare in solid state systems. Jiao *et al.* (Vidya Madhavan group, University of Illinois) performed STM measurements on the recently discovered heavy fermion superconductor  $\text{UTe}_2$  ( $T_c = 1.6$  K) and found clear signatures of competing Kondo effect and superconductivity within one unit cell. They also revealed ‘chiral’ in-gap states at step edges, which are theoretically predicted to exist at the boundaries of topological superconductors. Their result suggests that  $\text{UTe}_2$  is a strong candidate for chiral triplet topological superconductivity.

(a) STM image of  $\text{UTe}_2$ .  
(b,c)  $dI/dV$  spectra taken at the location denoted in (a), showing the Kondo resonance and the superconducting gap, respectively.  
(d) STM image of  $\text{UTe}_2$  with a terrace.  
(e)  $dI/dV$  spectra measured at the step edge marked by dots in (d), showing the chiral in-gap states.

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Ultra High Vacuum / Low Temperature / High Magnetic Field SPM System

Compact cryogenic high magnetic field SPM system

**Application fields:** monolayer films, topological matters, transition metal dichalcogenides etc.

## USM1500 Publication List in 2020 (Selected) 論文リスト

**"Packing Sierpiński Triangles Into Two-Dimensional Crystals"**

Y. Zhang *et al.*, J. Am. Chem. Soc. **142**, 17928 (2020).

**"Te-Vacancy Induced Surface Collapse and Reconstruction in Antiferromagnetic Topological Insulator MnBi<sub>2</sub>Te<sub>4</sub>"**

F. Hou *et al.*, ACS Nano **14**, 11262 (2020).

**"Kinetics-Limited Two-Step Growth of Van Der Waals Puckered Honeycomb Sb Monolayer"**

Z. Shi *et al.*, ACS Nano **14**, 16755 (2020).

**"Charge Density Modulation and the Luttinger Liquid State in MoSe<sub>2</sub> Mirror Twin Boundaries"**

Y. Xia *et al.*, ACS Nano DOI: 10.1021/acsnano.0c05397

**"Distinct Topological Surface States on the Two Terminations of MnBi<sub>4</sub>Te<sub>7</sub>"**

X. Wu *et al.*, Phys. Rev. X **10**, 031013 (2020).

**"Atomic-Scale Observation of Topological Vortices in the Incommensurate Charge Density Wave of 2H-TaSe<sub>2</sub>"**

S. Lim *et al.*, Nano Lett. **20**, 4801 (2020).

**"Tuning the Electronic Structure of an  $\alpha$ -Antimonene Monolayer Through Interface Engineering"**

Z. Shi *et al.*, Nano Lett. **20**, 8408 (2020).

**"Direct Observation of Full-Gap Superconductivity and Pseudogap in Two-Dimensional Fullerides"**

M. Ren *et al.*, Phys. Rev. Lett. **124**, 187001 (2020).

**"Reversible Potassium Intercalation in Blue Phosphorene–Au Network Driven by an Electric Field"**

Y. Liu *et al.*, J. Phys. Chem. Lett. **11**, 5584 (2020).


**"Heterostructures of Tellurium on NbSe<sub>2</sub> from Sub-Monolayer to Few-Layer Films"**

J. Xian *et al.*, Nanoscale **12**, 1994 (2020).

## "Packing Sierpiński Triangles into Two-Dimensional Crystals"

Y. Zhang *et al.*, J. Am. Chem. Soc. **142**, 17928 (2020). Product used: USM1500

Two-dimensional (2D) fractal patterns such as the Sierpiński triangles (STs) are theoretically predicted to exhibit intriguing optical and electronic properties due to their aperiodic but ordered structures. However, realizing large-scale 2D crystals has remained challenging. Zhang *et al.* (Yongfeng Wang group, Peking Univ.) succeeded in constructing a series of 2D packed crystals of STs consisting of 1,3-bis (4-pyridyl) benzene (BPyB) molecules and Fe atoms on Au surfaces. Their STM characterization combined with theoretical calculations suggests that molecular free diffusion and structure matching between STs and Au play important roles in the construction of 2D crystal of STs.

(a) The ST unit consisting of three BPyB and one Fe atom. (b) STM image of a typical ST on Au(111). (c-e) Large-scale STM image of 2D crystals of STs. The coverage is 1, 0.65 and 0.57 monolayer, respectively. (f-h) Close-up STM image of (c-e), respectively.

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Dilution-Refrigerator Based UHV/Ultra Low Temperature/High Magnetic Field STM System

Electron temperature ~110mK

**Application fields:** topological matters (Majorana), etc.

## USM1600 Publication List in 2020 (Selected) 論文リスト

**"Robust Hot Electron and Multiple Topological Insulator States in PtBi<sub>2</sub>"**

X. Nie *et al.*, ACS Nano **14**, 2366 (2020).

**"Localized Spin-Orbit Polaron in Magnetic Weyl Semimetal Co<sub>3</sub>Sn<sub>2</sub>S<sub>2</sub>"**

Y. Xing *et al.*, Nat. Commun. **11**, 5613 (2020).

**"A Tunable and Unidirectional One-Dimensional Electronic System Nb<sub>2n+1</sub>Si<sub>n</sub>Te<sub>4n+2</sub>"**

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Joule-Thomson Ultra High Vacuum Very Low Temperature High Magnetic Field SPM System

Low temperature/high magnetic field SPM with optical access

Technology taken over from SPECS. GmbH and improved and integrated with USM1200 series

**Application fields:** molecules, graphene, etc.

## JT-SPM Publication List in 2020 (Selected) 論文リスト

**"Resolving Quinoid Structure in Poly-Para-Phenylene Chains"**

B. Yuan *et al.*, J. Am. Chem. Soc. **142**, 10034 (2020).

**"Precise Control of  $\pi$ -Electron Magnetism in Metal-Free Porphyrins"**

Y. Zhao *et al.*, J. Am. Chem. Soc. **142**, 18532 (2020).

**"Atomically Precise Synthesis and Characterization of Heptaauthrene with Triplet Ground State"**

X. Su *et al.*, ACS Nano **20**, 6859 (2020).

**"Designer Spin Order in Diradical Nanographenes"**

Y. Zheng *et al.*, Nat. Commun. **11**, 6076 (2020).


**"Engineering of Magnetic Coupling in Nanographene"**

Y. Zheng *et al.*, Phys. Rev. Lett. **124**, 147206 (2020).

## "Precise Control of $\pi$ -Electron Magnetism in Metal-Free Porphyrins"

Y. Zhao, K. Jiang *et al.*, J. Am. Chem. Soc. **43**, 18532 (2020). Product used: JT-SPM

Porphyrin is one of the most fundamental building blocks for functional molecular systems. The magnetism of localized spins in the center of the porphyrin has long been studied, but realizing delocalized spins in the porphyrin, which offer potentials for building correlated quantum spins, has remained elusive. Zhao *et al.* (Wang, Zhuang, Jia group, Shanghai Jiao Tong Univ.) revealed that metal-free porphyrins adsorbed on Au(111) host delocalized singlet  $\pi$ -electron spins using STM/AFM and theoretical calculations. The delocalized magnetism was precisely controlled by tailoring peripheral  $\pi$ -electron topology in

metal-free porphyrins and reversely switched on/off via STM manipulation by tuning the interfacial charge transfer between the porphyrin and Au substrate. Their study provides an effective way to engineer  $\pi$ -magnetism induced by porphyrins, implying potential applications for organic-molecules-based quantum information and spintronics.

Magnetism switching experiment by STM manipulation. (a-c) From left to right: chemical structure, nc-AFM image with a CO tip, constant-height STM image. (d)  $dI/dV$  spectra taken at the marked location in STM images.


## Original New Functions

- ✓ **TCSPC fluorescence lifetime**
- ✓ Digital response capability **New**
- ✓ One-stroke scanner for thin film **New**
- ✓ Pseudo-CW performance enhanced

TCSPC option

Scanner software is incorporated in picoTAS

Digital response makes it easier to detect signal!

Demo Room for Optical Instruments Launched!

We Welcome Your Visit for Experiments

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
Conventional UV/VIS/NIR Flash Photolysis System

Normal Mode High Sens. Mode

Significant improvement of sensitivity!

- ✓ Digital response capability **New**
- ✓ Complete removal of Q-SW noise by pump-shutter **New**
- ✓ High sensitivity achieved by 5 Hz shutter, 12 bit oscilloscope and chopping mode **New**

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Cryostat for Spectrophotometer USP-203 Series

- ✓ Total number of shipments of CoolSpeK since released > 450
- ✓ Wide input voltage available now. No transformer required **New**
- ✓ Newly compatible with Cary 100 / 300 / Eclipse (Agilent Tech.), Spectrum Two (PerkinElmer) UH5700 (Hitachi)



Upcycled HP8453 by OLIS Inc. for CoolSpeK is now available from us.

## "Effect of Deuteration on Relaxation Dynamics of the Perylene Excimer Studied by Subnanosecond Transient Absorption Spectroscopy"

Y. Shibasaki, R. Katoh et al., J. Phys. Chem. A 125, 1359 (2021). Product used: picoTAS

Prof. Katoh's group in Nihon University found that the deuterated perylene excimer in solution had a longer lifetime than the undeuterated excimer. They studied the effect of deuteration on the relaxation dynamics of the excimer of perylene in solution using subnanosecond time-resolved transient absorption spectroscopy based on the randomly interleaved pulse-train (RIPT) method, and they concluded that the rate of internal conversion was suppressed by deuteration. Such a fundamental finding has been overlooked because it had been difficult to perform systematic study for phenomena lying in nanosecond region by conventional techniques. Prof. Katoh says, "This is actually a modest job. However, I believe that this is an important data and it would appear in textbooks in the future."

# Ying-Shuang Fu

Huazhong University  
of Science and Technology, China

## Research Interests

- Molecular beam epitaxy growth of low dimensional materials
- Electronic properties of topological states of matter
- Electronic properties of correlated states in low dimensional systems
- Spin-resolved spectroscopic imaging of magnetism in low dimensions


USM1300

0.4K-12T STM

USM1500


2K-8T STM/AFM

USM1200


5K STM

## Phason-Polaron Spectrum of Pure One-Dimensional Charge Density Wave in a Single $\text{Mo}_6\text{Se}_6$ Nanowire


X. Yang, *et al.* Phys. Rev. X, **10**, 031061 (2020)

## Spin Resolved Intrinsic Antiferromagnetic Order in a Van Der Waals Monolayer Crystal


## Topological Edge State Residing at Step Edge of $\text{WTe}_2$

L. Peng, *et al.*, Nat. Commun. **8**, 659 (2017).

## Electronic States of Bi Ultrathin Films Modulated by a $\text{NbSe}_2$ Substrate

L. Peng, *et al.*, ACS Nano **13**, 1885 (2019).

# Peter Liljeroth

Aalto University, Finland


## Research Interests

- Low-temperature scanning tunneling microscopy and spectroscopy
- Atomic force microscopy
- Two-dimensional materials, metal-organic frameworks
- On-surface synthesis
- Artificial lattices
- Designer quantum materials

## Artificial Atomic Lattices


### Examples of Artificial Lattices Constructed Using Chlorine Vacancies on $\text{Cu}(100)$ .

(a) STM topography, LDOS map and spatially resolved  $dI/dV$  spectroscopy on a dimer chain with two domain walls hosting topological mid-gap states. (b) STM topography (top) and the dispersive and flat band LDOS maps of a 1D cross chain. (c) Topography and LDOS maps of a Lieb lattice structure assembled from Cl vacancies.


**Topological States in Engineered Atomic Lattices**, R. Drost *et al.* Nat. Phys. **13**, 668 (2017).  
**Tunable Topological Domain Wall States in Engineered Atomic Chains**, Md N. Huda *et al.* npj Quant. Mater. **5**, 17 (2020).  
**Designer Flat Bands in Quasi-One-Dimensional Atomic Lattices**, Md N. Huda *et al.* Phys. Rev. Res. **2**, 043426 (2020).

## Yu-Shiba-Rusinov States and Spin-Flip Excitations in Molecular Systems

## Topological Superconductivity in Van Der Waals Heterostructures

### Topological Superconductivity in $\text{CrBr}_3$ - $\text{NbSe}_2$ Heterostructures.

(a) Schematic of the heterostructure. (b) STM topography and  $dI/dV$  spectroscopy showing the presence of edge modes at the edges of the  $\text{CrBr}_3$  monolayer. (c) Spatially resolved spectroscopy and LDOS maps of the Majorana zero modes.

**Topological Superconductivity in a Van Der Waals Heterostructure**, S. Kezilebieke *et al.* Nature **588**, 424 (2020).  
**Moiré-Enabled Topological Superconductivity**, S. Kezilebieke *et al.* arXiv:2011.09760.

(a),(b) STM topography and  $dI/dV$  spectroscopy of MnPc and CoPc molecules. (c) Tunneling conductance as a function of the exchange coupling  $pJ$  showing the co-existence of YSR states and spin-flip excitations.

**Coupled Yu-Shiba-Rusinov States in Molecular Dimers on  $\text{NbSe}_2$** , S. Kezilebieke *et al.* Nano Lett. **18**, 2311 (2018).  
**Observation of Coexistence of YSR States and Spin-Flip Excitations**, S. Kezilebieke *et al.* Nano Lett. **19**, 4614 (2019).

# Toru Hirahara

Tokyo Institute of Technology,  
Japan

## Research Interests

- Novel atomically thin 2D systems
- High-temperature superconductivity: FeSe/STO
- Magnetic topological insulators

# Yi Du

Australian Institute for Innovative  
Materials (AIIM)  
The University of Wollongong,  
Australia

## Research Interests

- Scanning tunneling microscopy
- 2D Xenes materials (silicene, germanene, stanene and blue phosphorene) and 2D frustrated lattice systems
- MBE growth of 2D device structures (photodetectors, sodium ion batteries, etc.)


## Interfacial Superconductivity in FeSe Ultrathin Films on SrTiO<sub>3</sub> Probed by *In Situ* Independently Driven Four-Point-Probe Measurements

A. Pedersen *et al.*, Phys. Rev. Lett. **124**, 227002 (2020).

### USM1400

Photograph of the USM1400 system

- (a) Photograph of the four probes. (b) Probe spacing dependence measurements to determine the conductivity dimension. (c) Temperature dependence of the measured resistance for films with different thicknesses. (d) Same as (c) but shown in normalized resistance to show the universal feature of the superconductivity onset at 40 K.

## Superconductivity of Single Unit Cell FeSe/SrTiO<sub>3</sub>(001): Substrate-Surface Superstructure Dependence

T. Tanaka *et al.*, Phys. Rev. B **98**, 121410(R) (2018).


## Superconducting Dome Revealed by Surface Structure Dependence in Single Unit Cell FeSe on SrTiO<sub>3</sub>(001)

T. Tanaka *et al.*, Phys. Rev. B **101**, 205421 (2020).

### USM1500

Photograph of the USM1500 system with the load lock, sample preparation, and measurement chambers


- (a)-(d) Atomically resolved image of the 1UC FeSe grown on different STO surface superstructures. (e):  $dI/dV$  spectra taken on the 1UC FeSe grown on different STO surface superstructures. (f) Same as (e) but those obtained near the Fermi level showing the superconducting gap.

## STM Facilities in the Team


## Realization of Flat Band with Possible Nontrivial Topology in Electronic Kagome Lattice

- (a) STM image of the Kagome lattice on a twisted multilayered silicene on the Ag(111) substrate. (b,c) STS image simultaneously obtained with (a) at the flat-band energy, showing the Kagome pattern and at the edge state energy, respectively.

Z. Li *et al.*, Sci. Adv. **4**, eaau4511 (2018).


## Experimental Realization of Two-Dimensional Buckled Lieb Lattice

H. Feng, *et al.*, Nano Lett. **20**, 2537 (2020).

- (a) Atomic resolution STM image of 0.5 ML Sn on Al(100). (b) Model structure of buckled Lieb lattice comprising Sn in a  $\sqrt{2} \times \sqrt{2}$  superstructure and the topmost Al atoms of the Al(100) substrate. (c) 3D band structure of the buckled Lieb lattice, possessing a broken nodal line loop.

## Realization of Strained Stanene by Interface Engineering

- (a) STM image of 1 ML Sn-deposited surface on the Au(111) substrate, showing a  $\sqrt{3} \times \sqrt{7}$  superstructure. (b,c)  $dI/dV$  spectra and in-situ Raman spectra for the  $2 \times 2$ ,  $\sqrt{3} \times \sqrt{3}$ , and  $\sqrt{3} \times \sqrt{7}$  phases, respectively.

Y. Liu *et al.*, J. Phys. Chem. Lett. **10**, 7, 1558 (2019).

The place where the first Japanese beer was brewed in Dojima

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### <Statistics>

Production of beer in Osaka is 2nd largest in Japan.  
Beer consumption per person in Osaka is 2nd largest in Japan.  
大阪のビールの生産量は国内2位、ビール一人当たり消費量2位

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### <History>

The grapes recently produced in Osaka originate from the transplant of Koshu grapes from Osaka nursery garden to Kashiwara city in 1878. Since then, the grape production was expanded to wide areas in the southern part of Osaka during the Meiji and Taisho Eras and thus has a long history of more than 100 years.

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The oldest winery in West Japan, founded by Sakujiro Takai in 1914  
The wine was served at G20 Osaka Summit 2019.  
高井作次郎創業(西日本最古) 2019年G20大阪サミットで提供

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# Where to Contact in Your Area 代理店紹介

■ SPM, ■ Optics, ■ SPM & Optics

UNISOKU has  
great distributors  
to take care of  
your inquiries in  
your area.

## EUROPE ISRAEL NORTH AFRICA

### nanoscore GmbH

Maisebachstraße 3, 61479  
Glashütten, Germany  
TEL +49 (0) 6174-619-950  
FAX +49 (0) 6174-619-949  
t.berghaus@nanoscore.de  
www.nanoscore.de

## RUSSIA

### LLC SPECS-TII RUS

5, Bolshaya Sadovaya Str.,  
Moscow, 123 001  
TEL +7 (495) 920-90-05  
info@specs-tii.ru  
www.specs-tii.ru

## CANADA USA MEXICO

### SPECS-TII Inc.

20 Cabot Blvd., Suite 300  
Mansfield, MA 02048 USA  
TEL +1-508-618-1292  
usa@specs-tii.com  
www.specs-tii.com

### Worldwide Exchange LLC

P. O. Box 1559  
Princeton, NJ 08542  
info@wwe-us.com  
www.wwe-us.com

## CHINA

### SPECS-TII Technology (Beijing) Co., Ltd.

Room 418A, Building 39, 4  
Workers Stadium North Rd.  
Chaoyang District, Beijing 100027,  
P.R.China  
TEL +86-10-65010355  
sales@specs-tii.com.cn  
www.specs-tii.com.cn

## TAIWAN

### SHUMOTEK CORP.

No. 111, Jiabei 1st St. Zhunan,  
Miaoli County 35058, Taiwan  
TEL +886-37-612810 (Office),  
+886-919-706759  
febchangs@gmail.com  
info@shumotek.com  
www.shumotek.com

## SOUTH KOREA

### INA Korea Co., Ltd.

12F, 110 Mallijae-ro,  
Mapo-gu, Seoul,  
04184 Korea  
TEL +82-2-714-5456  
AX +82-2-714-5457  
shpark@inacr.com  
www.inacr.com

### Lambda Ray Co., Ltd.

55-16, 247beon-gil,  
Mabuk-ro, Giheung-gu,  
Yongin-si, Gyeonggi-do,  
16891 Korea  
TEL +82-31-693-5776  
FAX +82-31-284-6089  
jayjeong@lambdaray.co.kr  
www.lambdaray.co.kr

For customers outside those areas, contact UNISOKU directly. Mail: info@unisoku.co.jp

## Events

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unisoku.co.jp

unisoku.com

UNISOKU Co., Ltd.

Web site: <http://www.unisoku.com/>

E-mail: [info@unisoku.co.jp](mailto:info@unisoku.co.jp)

2-4-3 Kasugano, Hirakata, Osaka 573-0131 Japan

TEL +81-72(858)6456 FAX +81-72(859)5655

20210510

## Global cooperation of SPECS-TII