2024 UNISOKU NEWSLETTER

UNISOKU Co., Ltd.

Web site: <https://www.unisoku.com/> E-mail: info@unisoku.co.jp

UNISOKU contributes to the development of science and technology by providing customers with measurement systems that meet their exploration needs.

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We have been constantly developing and manufacturing systems with a spirit of challenge to meet the needs of our customers. In the world of rapid technological progress, we will continue to strive to meet our customers' demands and to propose the best possible products.

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Hiroyuki Mizuno

We would like to thank our customers and all those involved in our business for the 50th year since our foundation. We aim to contribute to society and achieve steady growth by listening to our customers' needs and incorporating them into our management policy. We take our customers' opinions seriously and will continue to take on the challenge of realizing new 'uniqueness'.

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Yutaka Miyatake

The last three years have seen a series of global events that have forced everyone to change their previous values. We will continue to face our products and customers with a sense of gratitude and humbleness while considering what we can change and how we can contribute to society in the future.

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2023 Yearly Events

January 1

Time-resolved STM paper was published in Scientific Reports. UNISOKU and Univ. of Tsukuba conducted a joint press release.
時間分解STM論文がScientific Reportsに掲載。筑波大と共同でプレスリリース

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
A new member joined the Production department.
生産技術グループに新入社員入社

May 5

Development of USM1800 was featured in the Nikkan Kogyo Simbun.
USM1800開発が日刊工業新聞に掲載

The UNISOKU main building was refurbished and the dining room was newly set up.
本社を改装し食堂コーナーを設置

June 6

Three papers using UNISOKU systems were simultaneously published in Nature.
ユニソク装置を使用した論文が3本同時にNature掲載

UNISOKU contributed to ICP2023 in Sapporo as an exhibitor, an organizer of the luncheon seminar, and an oral presenter.
第31回光化学国際会議(札幌)において機器展示ランチョンセミナー口頭発表を実施

Founder Toshihiko Nagamura passed away.
創設者、長村俊彦氏が逝去


October 10

The Sales Meeting 2023 was held in UNISOKU.
セールスミーティング2023開催
A new member joined the Optical instruments department.
分光課に新入社員入社

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December 12

UNISOKU was awarded for The 57th (Fiscal Year 2023) Good Company Award Special Prize.

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History of our founder, Toshihiko Nagamura

Our founder, Toshihiko Nagamura, who retired in 2013, passed away on June 17, 2023 (age 86).
Once again, we would like to express our deepest gratitude to his foundation and great contribution to the development of the business, training of many employees, and field of science.

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1st UNISOKU STM USM-101
ユニソク最初のSTM

1st UNISOKU STM data
ユニソク最初のSTM像


Snapshot during a company trip to Korea
社員旅行、韓国にて

1989


Productized Ultra-high Vacuum STM systems and started its sales.
超高真空 STM 装置を製品化、販売開始

Ultra-high Vacuum STM USM-301U
超高真空 STM

1991

Productized Ultra-high Vacuum Low Temperature STM systems and started its sales.
超高真空低温 STM を製品化、販売開始

Ultra-high Vacuum Low Temperature STM system USM-501
超高真空低温 STM

Ultra-high Vacuum Low Temperature STM system USM-602
超高真空低温 STM

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1998

Built the new head office in Kasugano, Hirakata city (the current location) and moved head office.
枚方市春日野に新社屋が完成し、本社を移転

1999

Started sales of 2K High Magnetic Field STM, USM1300.
2K 磁場中 STM (USM1300初号機)の販売開始

2000

In the 2000s, the domestic SPM sales gradually dropped, at the same time, UNISOKU began to expand business in overseas.
2000年代に入るとSPM国内売上は緩やかに下落し、それと同時に海外展開を始めた。

2002

Productized 400mK Very Low Temperature USM1300 systems and started its sales.
400mK 極低温USM1300の製品化、販売開始

Delivered USM1200 to South Korea as the first STM in oversea.
海外向として初めてのSTM、USM1200を韓国に納入

Taken at a conference in China in 2006, he has taken the lead in developing SPM's overseas business.
2006年中国の学会にて、自ら先頭に立ちSPMの海外事業展開に注力

Awarded for "Technological Achievement Award" from The Japan Society for Analytical Chemistry for the development of fast reaction analysis/nanoscale surface analysis.

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
2003

Delivered USM1300 to Prof. Goran Lab of Argonne National Laboratory as 1st STM in the USA and joined its installation in person.

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2023 was really a sad year that we lost Nagamura-san.
I was very fortunate to meet Nagamura-san during my first trip to Japan.
He was a remarkable person with a big heart.
I cherish and remember all the moments that I spent with him in Japan and in Chicago.

From Prof. Goran Karapetrov, Drexel University

At RHK Sales Meeting in 2004
2004年RHK Sales Meetingにて

2004

Started sales of Ultra-high Vacuum Low Temperature SPM systems (USM-1400 series).

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Delivred USM1300 to Prof. Xue Lab of Chinese Academy of Sciences as 1st STM in China.

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Delivered 1st USM1400 in China with optical detector to Prof. Dong Lab of USTC.

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I am very sad to hear that Nagamura-san, my dear old friend, passed away two days ago. Yes, he has been a very good friend of mine all these years, with a lot of beautiful memories in developing advanced photon-STM systems together, plus happy times in enjoying life and nature. Please send my condolences to his family members.

May Nagamura-san rest in peace in heaven! He will live in our heart forever!
Sincerely,

From Prof. Zhen-Chao Dong, USTC

2007

Awarded Osaka Persons of Merit for New Technology Development Award from Osaka Prefecture for achievements in the invention of Lower Temperature Scanning Probe Microscope systems.

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2010

UNISOKU joined TII group and he was appointed as Chairman.

Shoji Suruga was iappointed as the second President.

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With the second President Suruga
二代目社長駿河と

At Team Appreciation Party in 2010
2010年慰労会にて

Awarded at the Surface Science Society of Japan for his work of Ultra High Vacuum Low Temperature Scanning Probe Microscope.
日本表面科学会 第一回産業賞 受賞

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1st USM1600 System with 40mK specification, which he greatly contributed to its depevelopment, was delivered to Prof. Xue's lab of Tsinghua University, China.

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His visit to the 2nd factory of UNISOIKU
新工場見学

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Referred from 45th Anniversary History of UNISOKU, News Letter 2019

"Since UNISOKU was founded, our motto has been 'Contribute to the progress in science and technology by developing unique instruments in collaboration with customers'. I hope UNISOKU will continue this policy to respond to expectations from customers all around the world."

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(Founder and Former Chairman) Toshihiko Nagamura

Sales Meeting 2023

Inviting nine people of seven distributors from the USA, China, Germany, Taiwan and Korea this year, UNISOKU finally held sales meeting for the first in four years!
As this event was originally scheduled to be held every two years, there was no end to the topics of conversation with the distributors about recent developments, such as business inquiries.

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This time, we were reminded of the importance of face-to-face meetings in which we could see facial expressions and reactions of the participants. After the presentation from UNISOKU, there were many questions. We are happy to see that our information is being conveyed more accurately to the distributors. We are sure that they will take what they learned from this meeting and apply it to their customer service in their own countries.

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
On the last day, we had an optional sightseeing day, which included a cruise around Osaka City on an amphibious bus and a winery tour at Kawachi Winery, which was featured in Newsletter 2021.

For even Osaka local people, it is rare to have a chance to ride on an amphibious bus. The view of Osaka city from the amphibious bus was refreshing. Kawachi Wine was once served at the G20 Osaka Summit banquet. The winery has been cultivating grapes since the mid-Meiji period and its taste of wine is characterized by its sweet aroma and dry taste.

At the winery, about 10 kinds of wine and plum wine were served for tasting. The wines were so good that participants bought them as souvenirs.

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Zum Greifenstein 5, 65594
Runkel, Germany
sales@nanoscore.de
www.nanoscore.de

Introduction of Dr. Andreas Bettac

We are proud to announce that Dr. Andreas Bettac has joined nanoscore. Andy substantially strengthens our sales and customer service activities. He brings in his decades-long experience in SPM instrumentation and market relations. In parallel with the growing number of products and applications this will greatly improve our market presence – and our ability of timely answers to customers' request.

SOUTH KOREA

INA Korea Co., Ltd.

12F, 110 Mallijae-ro,
Mapo-gu, Seoul, 04184
South Korea
shpark@inacr.com
www.inacr.com

Lambda Ray Co., Ltd.

23-13, Mabuk-dong,
Gihung-gu, Yongin-si,
Gyeonggi-do, 16891,
South Korea
jayjeong@lambdaray.co.kr
www.lambdaray.co.kr

CHINA

SPECS-TII Technology (Beijing) Co., Ltd.

Room 418A, Building 39,
4 Workers Stadium North Rd,
Chaoyang District, Beijing
100027, P.R.China
sales@specs-tii.com.cn
www.specs-tii.com.cn

TAIWAN

SHUMOTEK CORP.

No. 111, Jiabei 1st St. Zhunan,
Miaoli County 35058, Taiwan
febchangs@gmail.com
info@shumotek.com
shumotek.com.tw/

● SPM ● Optics ● SPM & Optics


UNISOKU is a member of TII Group, which is headed by our parent company, Tokyo Instruments, Inc. While respecting the core competencies of TII Group, we aim to "create new value" through strong cooperation.

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Redesigned most of optical system including monochromators in picoTAS Optics/Mechanics Unit (picoTAS-OMU)

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
- Redesign of optical system: improvement of chromatic aberrations, etc.
- Renewal of motor drive system: faster motion of most of motorized units
- Change of top cover: plastic to metal for better exogenous noise immunity (EMS)
- Adopted more compact pump laser in -ns model

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
-ns Model Lasers

Renewed Monochromator

Thermodynamic Control of Intramolecular Singlet Fission and Exciton Transport in Linear Tetracene Oligomers

S. Nakamura et al., Angew. Chem. Int. Ed. 62 e202217704 (2023).


Publication Introduction 論文紹介

Recently, singlet fission (SF), a multiexciton generation process that can raise the efficiency of exciton generation up to 200 %, has been actively investigated as a next-generation photo-energy utilization. To enhance the efficiency of SF,

Dr. Nakamura of Prof. Hasobe's group, Keio University and co-researchers proposed a molecular design strategy by linear oligomers, and they newly synthesized a series of homo-tetracene [(Tc)_n], and hetero-tetracene oligomers [TcF₃-(Tc)₄-TcF₃], then evaluated those performances by electron paramagnetic resonance (EPR) and femtosecond to millisecond transient absorption (TA) measurements. Especially, they made full use of CoolSpeK with picoTAS and TSP-2000, namely, they measured temperature dependence of rate constants of each process by TA to precisely calculate thermodynamic parameters and succeeded in explaining the high efficiency of SF in both (Tc)_n and TcF₃-(Tc)₄-TcF₃. Especially, efficient localization of exciton trapping was confirmed in TcF₃-(Tc)₄-TcF₃, with a trapped triplet yield of 176% due to the increase of entropy change, ΔS. Such thermodynamic control of SF by this molecular design will provide a new perspective for novel photofunctional systems such as quantum information science and solar energy conversion.

Illustration and Data Courtesy from Prof. Hasobe

(a) Chemical structures of Tc-oligomers, (b) Conceptual schemes, (c) psTA spectra of (Tc)₂ in toluene, (d) psTA spectra of TcF₃-(Tc)₄-TcF₃ in PhCN, (e) TA Decay profiles of (Tc)₄ at 530 nm at different temperatures, (f) The plots of k vs. temperature in TcF₃-(Tc)₄-TcF₃

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The 42nd Issue of 'Essentials in Chemistry' 化学の要点シリーズ42

In "Time-resolved measurement of chemical reactions with pulsed lasers: Transient Absorption Measurement", the 42nd issue of 'Essentials in Chemistry' (in Japanese), an explanatory article of RIPT method was written by us. 化学の要点シリーズ42『パルスレーザーによる化学反応の時間分解計測: 過渡吸収測定』の中で、RIPT法に関する解説記事執筆を担当。

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In the December 2023 issue of "KOUKAGAKU" (photochemistry)" (in Japanese), a magazine published from the Japanese Photochemistry Association, our director, Dr. Nakagawa, wrote the preface, "Supporting members that give opinions" 光化学協会の会誌2023年12月号にて、巻頭言執筆を担当。

CoolSpeK Updates CoolSpeKの最新情報

CoolLink Automatic Temperature Variable Software

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Automatic temperature-dependent spectra measurement is now available!

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
- Easy to design temperature profile with PC
- Monitorable actual temperature in real time
- Linkable with various commercial spectroscopy system

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Adapter for Agilent Cary 3500 Flexible Module

Agilent Cary3500 Flexible Moduleに対応


Adapter for Agilent Cary 3500 Flexible Module has been released. For more information, please have a look at the instruction movie on YouTube.

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CoolSpeK SLIM

CoolSpeK SLIM (USP-203C-ST-BP)

Cryostat for 2 mm light-path cuvette. It is suitable for picoTAS, pump-probe spectroscopy system, etc.

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CoolSpeK with Large Optical Windows

CoolSpeK 広角型 (USP-203C-SH-CD)

This model has larger optical windows than standard type has.* It is suitable for measuring circular dichroism spectrum. Most of accessories for standard type are also available with this model.

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We Now Offer In-House Experimental Demonstrations.

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Because we aim for after-purchase satisfaction, we provide our customers the opportunity to check the product performance before purchase. Further, we also aim to help customers who are not ready to purchase our systems conduct their research. To these ends, we have set up a room showing the newest instruments, both for demonstration purposes and for in-house experiment service.

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
Combined System of
Picosecond Transient Absorption and
TCSPC Fluorescence Lifetime

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We also offer on-site CoolSpeK demonstration at your facility. CoolSpeK adaptation to your spectrometer for custom demonstration measurements is available (only domestic)

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【Custom demo measurements】

We organize demonstration measurements of your samples using the HEMTO-TDS at our facility. Contact us to discuss the details of the samples you are interested in!

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
New version compatible with flag-type sample holders released!

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

STM image of HOPG

³He operation is now available!


- Under ⁴He operation, the base temperature is 1.1 K.
- Under ³He operation, the electron temperature is 0.9 K.

Time-Resolved Atomic Force Microscopy

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Our compact optical pump-probe system can be applied to AFM!

1. Tuning fork type time-resolved frequency modulation AFM


H. Mogi *et al.*, Appl. Phys. Express 17, 015003 (2024).

Collaboration with Prof. Shigekawa group
(Univ. of Tsukuba)

The introductory article was published.
“Next-Generation Time-Resolved Scanning
Probe Microscopy”
K. Iwaya *et al.*, Imaging & Microscopy 26, 34
(2024).

2. Time-resolved AFM in ambient conditions using a versatile SPM system

Collaboration with Dr. Minato
(Institute for Molecular Science)

Time-resolved force signals were
clearly detected!

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We are pleased to announce the launch of our 'Rental Lab' service, providing the necessary measurement environment for Ultra-low Temperature SPM measurements. To facilitate widespread and innovative use of Low-temperature SPM, we have established an ultra-low temperature SPM measurement platform.

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We have installed our high-end Ultra-low-temperature SPM equipment within our company premises and offer measurement environments to researchers worldwide on a paid basis. Our services include on-site experiments and remote experiments via online connections, tailored to meet your specific requirements.

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

Available Now

Carrier dynamics measurement of micro samples on insulating substrate

Demo measurement conditions

- Temperature: 77 K or 300 K
- Pressure: $\sim 10^{-6}$ Pa
- Laser wavelength: 488, 532 nm
- Temporal resolution: ~ 80 ps (532 nm), ~ 10 ns (488 nm)

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USM1600

Specifications

- $T_{STM\ Head} = 40mK$
- Vector Magnet operation
- RF-STM
- Long-term dl/dV measurement

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USM1200 JT

Specifications

- $T_{STM\ Head} = 1.5$ K (when optical shutters close)
- Compatible with AFM measurement
- Optical access capabilities by inertial-driven lens stages
- Time resolved STM with high spatial resolution
- Shot noise measurement by integrated RydeenAmp

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
Introduction of Publications

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Publication Stats in 2023

- Total number of publications using UNISOKU systems = 314 (294 in 2022)
- Total impact factors ~ 2610 (2268 in 2022)
Corresponding to 40 Nature papers (45 in 2022)
c.f. impact factor of Nature ~ 64 (50 in 2022)
- Impact factor per employee ~ 51 (~ 47 in 2022)

Popular Research Fields	Num. of Publications	Average Impact Factor
Low Dimensional Materials excluding TMDs, graphene, 2D superconductivity	36	8.8
Molecules	29	11.5
Transition Metal Dichalcogenides (TMDs)	27	14.2
Graphene (Twisted Bilayer Graphene)	22	14.2
Kagome Materials	20	14.0
2D Superconductivity	18	9.6
Topological Materials (Majorana, Weyl)	14	11.2
Fe-based Superconductors	14	13.5
Single Atom Spin (ESR-STM)	7	23.2

Publication List in 2023

Nature

1. Smectic Pair-Density-Wave Order in $\text{EuRbFe}_2\text{As}_4$
H. Zhao *et al.*, Nature **618**, 940 (2023). USM1300
2. Magnetic-Field-Sensitive Charge Density Waves in the Superconductor UTe_2
A. Aishwarya *et al.*, Nature **618**, 928 (2023). USM1300
3. Detection of a Pair Density Wave State in UTe_2
Q. Gu *et al.*, Nature **618**, 921 (2023). USM1500
4. Imaging Inter-Valley Coherent Order in Magic-Angle Twisted Trilayer Graphene
H. Kim *et al.*, Nature **623**, 942 (2023). USM1300

Science

- An Atomic-Scale Multi-Qubit Platform
Y. Wang *et al.*, Science **382**, 87 (2023). USM1300

Nature Materials

1. Tip-Induced Excitonic Luminescence Nanoscopy of An Atomically-Resolved Van Der Waals Heterostructure
L. Parra Lopez *et al.*, Nat. Mater. **22**, 482 (2023). USM1400
2. Tunable Spin and Valley Excitations of Correlated Insulators in Γ -Valley Moiré Bands
B. Foutty *et al.*, Nat. Mater. **22**, 731 (2023). USM1300

Advanced Materials

- Spin-Resolved Imaging of Antiferromagnetic Order in Fe₄Se₅ Ultrathin Films on SrTiO₃**
W. Zhang *et al.*, Adv. Mater. **35**, 2209931 (2023). [USM1500](#)
- Evidence of Nodal Superconductivity in Monolayer 1H-TaS₂ with Hidden Order Fluctuations**
V. Vano *et al.*, Adv. Mater. **35**, 2305409 (2023). [USM1300](#)

Advanced Energy Materials

Origin of Enhanced Overall Water Splitting Efficiency in Aluminum-Doped SrTiO₃ Photocatalyst
D. Murthy *et al.*, Adv. Energy Mater. **13**, 2302064 (2023). [USP-0BB-75XE](#)

Chem

Photodynamic Treatment of Acute Vascular Occlusion by Using an Iron–Nitrosyl Complex
J. Choe *et al.*, Chem **9**, 1309 (2023). [CoolSpeK](#)

Nature Chemistry

- Air- And Photo-stable Luminescent Carbodicarbene-Azaboracene Ions**
C. Deng *et al.*, Nat. Chem. **16**, 437 (2023). [CoolSpeK](#)
- Quantum Nanomagnets in On-Surface Metal-Free Porphyrin Chains**
Y. Zhao *et al.*, Nat. Chem. **15**, 53 (2023). [JT](#)

Nature Physics

- Evidence for Chiral Superconductivity On a Silicon Surface**
F. Ming *et al.*, Nat. Phys. **19**, 500 (2023). [USM1300](#)

Imaging Inter-Valley Coherent Order in Magic-Angle Twisted Trilayer Graphene

Kim, Choi *et al.*, Nature **623**, 942 (2023).

Product used: USM1300

This study focuses on magic-angle twisted trilayer graphene (MATTG), exploring its correlated electronic phases that break underlying symmetries. Using scanning tunneling microscopy, Kim, Choi *et al.* (Nadj-Perge group, California Institute of Technology) identify interaction-driven spatial symmetry breaking in low-strain samples, revealing atomic-scale reconstruction and correlated gaps in the tunneling spectrum. The observed Kekulé supercell suggests spontaneous inter-valley coherence between electrons, persisting across magnetic fields and temperatures associated with gap development. Large-scale maps show the coexistence of atomic-scale reconstruction and translation symmetry breaking at the longer moiré scale. The auto-correlation and Fourier analyses support the theoretically proposed incommensurate Kekulé spiral order. The findings provide insights into MATTG correlated phases, suggesting that superconductivity arises from an inter-valley coherent parent state.


Figure (a) dI/dV map at $V_{\text{Gate}} = 0 \text{ V}$, $+9 \text{ V}$ and its Fourier-transformed image showing V_{Gate} dependent Kekulé order on MATTG. **(b)** Evidence of moiré translation symmetry breaking. (Left) Large-scale dI/dV map including 20 moiré AAA sites. (Right) Kekulé order in Pattern A and B in the left image, showing inverted contrast between neighboring AAA sites.

Nature Physics-2

- Unidirectional Coherent Quasiparticles in the High-Temperature Rotational Symmetry Broken Phase of AV₃Sb₅ Kagome Superconductors**
H. Li *et al.*, Nat. Phys. **19**, 637 (2023). [USM1300](#)
- Unidirectional Electron-Phonon Coupling in the Nematic State in a Kagome Superconductor**
P. Wu *et al.*, Nat. Phys. **19**, 1143 (2023). [USM1300](#)
- Electronic Nematicity in the Absence of Charge Density Waves in a New Titanium-Based Kagome Metal**
H. Li *et al.*, Nat. Phys. **19**, 1591 (2023). [USM1300](#)
- Hofstadter States and Re-Entrant Charge Order in a Semiconductor Moiré Lattice**
C. Kometter *et al.*, Nat. Phys. **19**, 1861 (2023). [USM1300](#)

National Science Review

Structure Transformation from Sierpiński Triangles to Chains Assisted by Gas Molecules
C. Li *et al.*, Natl. Sci. **10**, nwad088 (2023). [USM1500](#)

Science Bulletin

- Charge Instability of Topological Fermi Arcs in Chiral Crystal CoSi**
Z. Rao *et al.*, Sci. Bull. **68**, 165 (2023). [USM1300](#)
- Inducing Itinerant Ferromagnetism by Manipulating van Hove Singularity in Epitaxial Monolayer 1T-VSe₂**
J. Zong *et al.*, Sci. Bull. **68**, 990 (2023). [USM1300](#)

ACS Nano-1


- Sensing the Local Magnetic Environment Through Optically Active Defects in a Layered Magnetic Semiconductor**
J. Klein *et al.*, ACS Nano **17**, 288 (2023). [USM1400-4P](#)
- The Bulk van der Waals Layered Magnet CrSBr is a Quasi-1D Material**
J. Klein *et al.*, ACS Nano **17**, 5316 (2023). [USM1400-4P](#)
- Single-Atomic-Layer Stanene on Ferromagnetic Co Nanoislands with Topological Band Structures**
C. Chen *et al.*, ACS Nano **17**, 7456 (2023). [USM1300](#)
- Realizing a Superconducting Square-Lattice Bismuth Monolayer**
E. Oh *et al.*, ACS Nano **17**, 7604 (2023).
- Inelastic Light Scattering in the Vicinity of a Single-Atom Quantum Point Contact in a Plasmonic Picocavity**
S. Liu *et al.*, ACS Nano **17**, 10172 (2023). [USM1400](#)
- Tuning Chirality of Self-Assembled PTCDA Molecules on a Au(111) Surface by Na Coordination**
Z. Liang *et al.*, ACS Nano **17**, 10938 (2023). [USM1400](#)

Evidence for Chiral Superconductivity on a Silicon Surface

Ming, Wu *et al.*, Nat. Phys. **19**, 500 (2023).

Product used: USM1300

Low-dimensional superconductivity has attracted considerable interest for decades. Tin (Sn) adatoms on a silicon (Si) substrate at a one-third monolayer coverage form a two-dimensional triangular lattice with unpaired electrons. The electrons organize into an antiferromagnetic Mott-insulating state, but doping the Sn layer with holes transforms it into a two-dimensional conductor that exhibits superconductivity at low temperatures. The repulsive interactions and frustration due to the triangular lattice suggest the possibility of unconventional superconductivity. Using scanning

tunneling techniques, Ming *et al.* (Wang group, Southern University of Science and Technology and Weitering group, Univ. of Tennessee) find a doping-dependent superconducting critical temperature, a fully gapped order parameter, evidence of time-reversal symmetry breaking, and enhanced zero-bias conductance near superconducting domain edges. These results hint at the possibility of Sn/Si(111) being an unconventional chiral d-wave superconductor.

Figure (a) STM image of $(\sqrt{3}\times\sqrt{3})$ -Sn surface on Si(111). **(b)** Doping dependence of the superconducting gap. $T = 0.5 \text{ K}$. **(c)** Conductance map at zero bias. **(d)** Fourier-transform of (c) showing a characteristic flower-like feature at $q = 0$. **(e)** dI/dV spectra across the substitutional Si defect showing two gap states (shaded).

ACS Nano-2

- 7. Double-Resonance Spectroscopy of Coupled Electron Spins on a Surface**
S. Phark *et al.*, ACS Nano **17**, 14144 (2023). [USM1300](#)
- 8. Influence of the Magnetic Tip on Heterodimers in Electron Spin Resonance Combined with Scanning Tunneling Microscopy**
X. Zhang *et al.*, ACS Nano **17**, 16935 (2023). [USM1300](#)
- 9. Optical Imaging of a Single Molecule with Subnanometer Resolution by Photoinduced Force Microscopy**
T. Yamamoto *et al.*, ACS Nano **18**, 1724 (2023). [USM1400](#)

Nature Communications-1

- 1. Charge Order Driven by Multiple-Q spin Fluctuations in Heavily Electron-Doped Iron Selenide Superconductors**
Z. Chen *et al.*, Nat. Commun. **14**, 2023 (2023).
- 2. Controllable Dimensionality Conversion between 1D and 2D CrCl₃ Magnetic Nanostructures**
S. Lu *et al.*, Nat. Commun. **14**, 2465 (2023). [USM1300](#)
- 3. Spectroscopic Signature of Obstructed Surface States in SrIn₂P₂**
X. Liu *et al.*, Nat. Commun. **14**, 2905 (2023). [USM1500](#)
- 5. Single-Electron Charge Transfer into Putative Majorana and Trivial Modes in Individual Vortices**
J. Ge *et al.*, Nat. Commun. **14**, 3341 (2023). [USM1500](#)
- 6. A Unique Van Hove Singularity in Kagome Superconductor CsV_{3-x}Ta_xSb₅ with Enhanced Superconductivity**
Y. Luo *et al.*, Nat. Commun. **14**, 3819 (2023). [USM1300](#)
- 7. Manipulating Single Excess Electrons in Monolayer Transition Metal Dihalide**
M. Cai *et al.*, Nat. Commun. **14**, 3691 (2023). [USM1500](#)
- 8. Towards Layer-Selective Quantum Spin Hall Channels in Weak Topological Insulator Bi₄Br₂I₂**
J. Zhong *et al.*, Nat. Commun. **14**, 4964 (2023). [SNOM1400](#)
- 9. From Stoner to Local Moment Magnetism in Atomically Thin Cr₂Te₃**
Y. Zhong *et al.*, Nat. Commun. **14**, 5340 (2023). [USM1300](#)

Single-Electron Charge Transfer into Putative Majorana and Trivial Modes in Individual Vortices

Ge *et al.*, Nat. Commun., **14**, 3341 (2023).

Product used: USM1500

This study investigates Majorana bound states, which are collective excitations in solids with properties of Majorana fermions, known for being their own antiparticles. In iron-based superconductors, zero-energy states in vortices have been suggested as potential Majorana bound states, but the evidence has remained debated. Using scanning tunneling noise spectroscopy, Ge *et al.* (Milan Allan group, Leiden University) examined tunneling into vortex bound states in the superconductors NbSe₂ and FeTe_{0.55}Se_{0.45}. They observed charge transfer of a single electron in both cases. For FeTe_{0.55}Se_{0.45}, the findings exclude Yu–Shiba–Rusinov states, supporting the possibility of both Majorana and trivial vortex bound states. The results provide insights into exotic states in vortex cores and implications for future Majorana devices, but further theoretical investigations are needed, especially regarding charge dynamics and superconducting tips.


Figure (a) Scanning tunneling noise microscopy setup. (b) Vortex lattice image of NbSe₂ and FeTe_{0.55}Se_{0.45}. (c) Local DOS across the vortex core in NbSe₂ (splitting) and FeTe_{0.55}Se_{0.45} (no-splitting). (d) Local noise spectroscopy on and off vortices in NbSe₂ and FeTe_{0.55}Se_{0.45}.

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Electric Control of Spin Transitions at the Atomic Scale

Kot *et al.*, Nat. Commun. **14**, 6612(2023).

Product used: USM1300

This study addresses the longstanding pursuit of electric spin control in solid-state physics, aiming to enhance information processing efficiency. The approach involves transitioning spintronics to the atomic scale. Using electron spin resonance scanning tunneling microscopy (ESR-STM), Kot *et al.* (Ast group, Max-Planck-Institut für Festkörperforschung) demonstrate electric control of spin resonance transitions in individual TiH molecules. Significant bias voltage-dependent shifts in the ESR signal, approximately ten times its line width, are observed. These shifts are attributed to the electric field in the tunnel junction, causing a displacement of the spin system and altering the g-factor and effective magnetic field of the tip. The study further demonstrates direct electric control of spin transitions in coupled TiH dimers, presenting new possibilities for coherent control of interconnected spin systems and contributing to an enhanced understanding of spin-electric coupling.


Figure (a) STM image of TiH molecules on 2ML MgO. (b) ESR signals at individual TiH molecule as a function of magnetic field and bias voltage. (c, d) Extracted g-factor and tip field as a function of bias voltage, respectively. (e, f) Interaction tuning in two different TiH molecule dimers.

- Collective Quantum Magnetism in Nitrogen-Doped Nanographenes**
G. Zhu *et al.*, *J. Am. Chem. Soc.* **145**, 7136 (2023). [JT-STM, AFM](#)
- Remote-Triggered Domino-like Cyclodehydrogenation in Second-Layer Topological Graphene Nanoribbons**
C. Ma *et al.*, *J. Am. Chem. Soc.* **145**, 10126 (2023). [USM1300](#)
- Revealing Intramolecular Isotope Effects with Chemical-Bond Precision**
X. Zhu *et al.*, *J. Am. Chem. Soc.* **145**, 13839 (2023). [USM1400TERS](#)
- Assembling Surface Molecular Sierpiński Triangle Fractals via K⁺-Invoked Electrostatic Interaction**
J. Dai *et al.*, *J. Am. Chem. Soc.* **145**, 13531 (2023). [USM1200](#)
- Realization of Long Operational Lifetimes in Vacuum-Deposited Organic Light-Emitting Devices Based on para-Substituted Pyridine Carbazolylgold(III) C⁺C⁺N Complexes**
C. Wong *et al.*, *J. Am. Chem. Soc.* **145**, 2638 (2023). [picoTAS](#)
- Reductive Coupling of Nitric Oxide by Cu(I): Stepwise Formation of Mono- and Dinitrosyl Species En Route to a Cupric Hyponitrite Intermediate**
M. Bhadra *et al.*, *J. Am. Chem. Soc.* **145**, 2230 (2023). [CoolSpeK](#)
- Generation, Spectroscopic Characterization, and Computational Analysis of a Six-Coordinate Cobalt(III)-Imidyl Complex with an Unusual S = 3/2 Ground State that Promotes N-Group and Hydrogen Atom-Transfer Reactions with Exogenous Substrates**
J. Yang *et al.*, *J. Am. Chem. Soc.* **145**, 26106 (2023). [CoolSpeK](#)
- Electronic Structure and Reactivity of Mononuclear Nonheme Iron–Peroxo Complexes as a Biomimetic Model of Rieske Oxygenases: Ring Size Effects of Macrocyclic Ligands**
W. Zhu *et al.*, *J. Am. Chem. Soc.* **146**, 250 (2023). [CoolSpeK](#)


Science Advances

Growth of Self-Integrated Atomic Quantum Wires and Junctions of a Mott Semiconductor
T. Asaba *et al.*, *Sci. Adv.* **9**, eabq5661 (2023). [USM1300](#)

Tunable Spin and Valley Excitations of Correlated Insulators in Γ -Valley Moiré Bands

Foutty *et al.*, *Nat. Mater.* **22**, 731 (2023). Product used: USM1300

This study investigates moiré superlattices formed by transition metal dichalcogenides, specifically focusing on twisted double-bilayer WSe₂. Using a scanning single electron transistor (SET), Foutty *et al.* (Feldman group, Stanford Univ.) conducted local electronic compressibility measurements to explore electronic correlations within the moiré bands associated with the Γ -valley. The study identifies charge-ordered phases at various integer and fractional moiré fillings. The magnetic field dependence of energy gaps and chemical potential measurements upon doping reveal spin-polarized ground states with spin-polaron quasiparticle excitations. Additionally, applying a displacement field

induces a metal-insulator transition by tuning between Γ - and K-valley moiré bands. The findings highlight the control over spin and valley characteristics of the correlated ground and excited states of twisted double-bilayer WSe₂

Figure
(a) Scanning SET setup.
(b) Inverse electronic compressibility $d\mu/dn$ as a function of filling factor ν and magnetic field B .
(c) Extracted charge gaps at $\nu = -1, -1/3$ as a function of B .
(d) Dual-gated capacitance device for displacement field tuning.
(e) $d\mu/dn$ as a function of filling factor ν and displacement field D .
(f) Measured gaps at $\nu = -1, -1/3$ as a function of D .

Small Fermi Pockets Intertwined with Charge Stripes and Pair Density Wave Order in a Kagome Superconductor
H. Li *et al.*, *Phys. Rev. X* **13**, 031030 (2023). [USM1300](#)

ACS Catalysis

- Intracluster O–O Coupling Pathway Evidenced for an Anderson-Type Single-Cobalt Polymolybdate Water Oxidation Catalyst**
N. Taira *et al.*, *ACS Catal.* **13**, 3211 (2023). [RSP-2000](#)
- Characterization and Reactivity of an Incredibly Reactive Intermediate in the Protonation Reaction of Dioxo-Manganese(V) Porphyrin with Acid**
Y. Katogi *et al.*, *ACS Catal.* **13**, 4842 (2023). [CoolSpeK](#)
- Reversible Deactivation of Manganese Catalysts in Alkene Oxidation and H₂O₂ Disproportionation**
J. Kasper *et al.*, *ACS Catal.* **13**, 6403 (2023). [CoolSpeK](#)
- Nonlinear Acid Promotion of Oxidation Reactions by Mononuclear Nonheme Iron(III)-Aqua Complexes**
M. Nilajakar *et al.*, *ACS Catal.* **14**, 34 (2024). [CoolSpek, RSP-601](#)

Proc. Natl. Acad. Sci. USA

First Order Quantum Phase Transition in the Hybrid Metal-Mott Insulator Transition Metal Dichalcogenide 4Hb-TaS₂
A. Nayak *et al.*, *PNAS* **120**, e2304274120 (2023). [USM1300](#)

Nano Letters -1

- Spatially Resolving Electron Spin Resonance of π -Radical in Single-molecule Magnet**
R. Kawaguchi *et al.*, *Nano Lett.* **23**, 213 (2023). [USM1300](#)
- Self-Intercalated 1T-FeSe₂ as an Effective Kagome Lattice**
Z. Zhang *et al.*, *Nano Lett.* **23**, 954 (2023). [USM1500](#)
- Oscillatory Order–Disorder Transition during Layer-by-Layer Growth of Indium Selenide**
Z. Chen *et al.*, *Nano Lett.* **23**, 1077 (2023). [USM1300](#)
- Real-Space Mapping of Local Subdegree Lattice Rotations in Low-Angle Twisted Bilayer Graphene**
Y. Ren *et al.*, *Nano Lett.* **23**, 1836 (2023). [USM1300](#)

Electronic Nematicity Without Charge Density Waves in Titanium-based Kagome Metal

Li *et al.*, *Nat. Phys.* **19**, 1591 (2023). Product used: USM1300

This study investigates layered crystalline materials with transition metal atoms arranged on a kagome network, focusing on the kagome metal CsTi₃Bi₅, which is isostructural to the vanadium-based kagome superconductor AV₃Sb₅. It is known that, unlike AV₃Sb₅, CsTi₃Bi₅ does not exhibit a detectable charge density wave state. Li *et al.* (Zeljko group, Boston College) performed spectroscopic-imaging scanning tunneling microscopy measurements and density functional theory calculations to reveal significant electronic correlation effects at low energies. They discover an electronic anisotropy breaking the sixfold symmetry of the lattice, originating from both in-plane and out-of-plane titanium-derived d orbitals. The findings suggest a hexagonal analog of nematic bond order in Fe-based superconductors, highlighting the role of electronic orbitals in CsTi₃Bi₅.


Figure
(a) STM image of Bi-terminated surface of CsTi₃Bi₅ showing a hexagonal lattice.
(b) Normalized conductance map showing a clear quasiparticle interference pattern.
(c) Fourier transform (FT) of (b) showing two-fold symmetry.
(d) Angle-dependent FT amplitudes of q_2 and q_3 in (c).
(e) Energy dispersion of q_2 and q_3 along A, B and C in (c).

Nano Letters -2

5. **Realization of Multiple Charge-Density Waves in NbTe₂ at the Monolayer Limit**
Y. Bai *et al.*, Nano Lett. **23**, 2107 (2023). [USM1300](#)
6. **Giant Periodic Pseudomagnetic Fields in Strained Kagome Magnet FeSn Epitaxial Films on SrTiO₃(111) Substrate**
H. Zhang *et al.*, Nano Lett. **23**, 2397 (2023).
7. **Observation of Robust and Long-Ranged Superperiodicity of Electronic Density Induced by Intervalley Scattering in Graphene/Transition Metal Dichalcogenide Heterostructures**
M. Zhang *et al.*, Nano Lett. **23**, 2630 (2023). [USM1300](#)
8. **Real-Space Observation of Unidirectional Charge Density Wave and Complex Structural Modulation in the Pnictide Superconductor Ba_{1-x}Sr_xNi₂As₂**
T. Qin *et al.*, Nano Lett. **23**, 2958 (2023). [USM1600](#)
9. **Visualization of Moiré Magnons in Monolayer Ferromagnet**
S. Ganguli *et al.*, Nano Lett. **23**, 3412 (2023). [USM1300](#)
10. **Control of the Magnetic Interaction between Single-Molecule Magnet TbPc₂ and Superconductor NbSe₂ Surface by an Intercalated Co Atom**
F. Ara *et al.*, Nano Lett. **23**, 6900 (2023). [USM1300](#)
11. **Kinetics of Nanobubbles in Tiny-Angle Twisted Bilayer Graphene**
C. Yan *et al.*, Nano Lett. **23**, 8532 (2023). [USM1500](#)
12. **Moiré Enhanced Two-Band Superconductivity in a MnTe/NbSe₂ Heterojunction**
J. Nie *et al.*, Nano Lett. **23**, 8370 (2023). [USM1300](#)
13. **Manipulating the Spin Orientation of Co Atoms Using Monatomic Cu Chains**
N. Noei *et al.*, Nano Lett. **23**, 8988 (2023). [USM1300](#)
14. **Quantum Phase Transition in Magnetic Nanographenes on a Lead Superconductor**
Y. Liu *et al.*, Nano Lett. **23**, 9704 (2023). [USM1300, JT](#)
15. **Chiral Charge Density Wave and Backscattering-Immune Orbital Texture in Monolayer 1T-TiTe₂**
M. Ren *et al.*, Nano Lett. **23**, 10081 (2023). [USM1500](#)
16. **Room-Temperature Ferromagnetism in Epitaxial Bilayer FeSb/SrTiO₃(001) Terminated with a Kagome Lattice**
H. Zhang *et al.*, Nano Lett. **24**, 122 (2023). [USM1300](#)
17. **Melting of Unidirectional Charge Density Waves across Twin Domain Boundaries in GdTe₃**
S. Lee *et al.*, Nano Lett. **23**, 11219 (2023). [USM1200](#)

Physical Review Letters

1. **Phase Shift and Magnetic Anisotropy Induced Field Splitting of Impurity States in (Li_{1-x}Fe_x)OHFeSe Superconductor**
T. Zhang *et al.*, Phys. Rev. Lett. **130**, 206001 (2023). [USM1600](#)
2. **Real-Space Imaging of Triplon Excitations in Engineered Quantum Magnets**
R. Drost *et al.*, Phys. Rev. Lett. **131**, 086701 (2023). [USM1300](#)
3. **Probing Hidden Mott Gap and Incommensurate Charge Modulation on the Polar Surfaces of PdCrO₂**
C. Wen *et al.*, Phys. Rev. Lett. **131**, 116501 (2023). [USM1200](#)
4. **Squeezed Abrikosov-Josephson Vortex in Atomic-Layer Pb Superconductors Formed on Vicinal Si(111) Substrates**
Y. Sato *et al.*, Phys. Rev. Lett. **130**, 106002 (2023). [USM1300](#)
5. **Oscillation of Electronic-Band-Gap Size Induced by Crystalline Symmetry Change in Ultrathin PbTe Films**
K. Chang *et al.*, Phys. Rev. Lett. **131**, 016202 (2023). [USM1600](#)
6. **Observation of Electronic Strong Correlation in VTe₂-2√3×2√3 Monolayer**
W. Zhao *et al.*, Phys. Rev. Lett. **131**, 086501 (2023). [USM1600](#)
7. **Iron Vacancy Tunable Superconductor-Insulator Transition in FeSe/SrTiO₃ Monolayer**
C. Xue *et al.*, Phys. Rev. Lett. **131**, 256002 (2023). [USM1600](#)

Chemistry A European Journal

- Generation and Characterization of a Tetraradical Embedded in a Curved Cyclic Paraphenylene Unit
Y. Miyazawa *et al.*, Chem. A Euro. J., **29**, e202301009, (2023). [picoTAS](#)

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Research Interests

- STM-induced Luminescence
- Tip-enhanced Photoluminescence and Raman Spectroscopy
- Tip-enhanced Photochemistry
- Single Molecules, Atomically-precise Graphene Nanoribbons and 2D Materials


USM1400-LT TERS


Features:

- Optimized Optical Access
- Time-Resolved Optical Detection (Hanbury Brown and Twiss Setup)
- Tunable-Laser Excitation

Tip-induced Excitonic Luminescence Nanoscopy of 2D Materials

Submolecular-Scale Control of Phototautomerization

Selected References:


- (1) L. Parra E. Lopez *et al.*, Nature Materials, **22** 482 (2023).
- (2) A. Roslowska *et al.*, Nature Nanotechnology, DOI: 10.1038/s41565-024-01622-4 (2024).

Kazuhiro Fujita and Abhay Pasupathy

Brookhaven National Laboratory, Upton NY, USA

Research Interests

- Spectroscopic Imaging Scanning Tunneling Microscopy
- High Temperature Superconductivity
- Electronic Interactions and Correlated Phases
- Topological States
- van der Waals Materials and Other 2D Systems
- Moiré Superlattices and Emergent Ground/Excited States


Shichao Yan

School of Physical Science and Technology, ShanghaiTech University, China

Research Interests

- Probing and Tuning Collective Electronic States in Layered Materials by Low-Temperature STM
- Atomic-Scale Spin Dynamics
- Development of Advanced STM Technique


USM1300


Features:

- Three Independent Contacts for Device Experiment
- 11 T Uniaxial Z-Magnet

Pair Density Wave

Spatially Modulated Gap in EuRbFe₄As₄

Selected References: (1) H. Zhao, *et al.*, Nature **618**, 940 (2023).


USM1600


Features:

- Ultra-Low Temperature ~30 mK
- 15 T Out-of-Plane Magnet

Unidirectional CDW and Complex Structure in NiAs-Based Superconductor

Incommensurate Charge Modulation in Delafossite Metals


USM1200


Features:

- Low-Temperature and Variable-Temperature STM
- Flexible for Laser-Coupled STM Research

Inducing and Tuning Kondo Screening in 1T-TaS₂

Water-Driven Reversible CDW Phase Transition in 1T-TaS₂

Selected References:

- (1) S. Shen *et al.*, Nature Communications **13**, 2156 (2022).
- (2) C. Wen, *et al.*, Physical Review Letters **126**, 256402 (2021).
- (3) C. Wen, *et al.*, Physical Review Letters **131**, 116501 (2023).
- (4) S. Shen, *et al.*, Nano Letters, **20**, 8854 (2020).
- (5) P. Kong, *et al.*, Nano Letters, **22**, 5635 (2022).
- (6) T. Qin, *et al.*, Nano Letters, **23**, 2958 (2023).

About Japanese Convenience Stores

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Convenience stores (Kombini) were first established in Japan in 1969. (There are various views.) Since 1975, they started 24-hour operation. Nowadays, the number of stores has increased to more than 57,000 across Japan and it is now common description to see "Convenience stores within 5 minutes walking distance" in real estate information.

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The floor area of the stores is very compact, and the products are arranged in typically the same layout in any stores. The selection is very extensive, offering lunchboxes, bread, sweets, beverages, alcohol, office supplies, and household goods. Especially for foods and beverages,

カップのお酒 内容量増量も同料金!!
Cup sake, upsized version is also available at the same price!

new products are introduced every week, and you may be able to find something that reflects the seasonal occasions. For example, in January many strawberry sweets were on display.

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You can also heat up the food you buy in the microwave oven installed in the store, and some stores have eat-in corners. When you have a sudden craving for snacks in the middle of the night or a cup of coffee made by an espresso machine in the early morning, that will come true.

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