

# Liquid Helium Level gage HLM-100

Displays the liquid level of Helium tank for USM 1400.  
Makes it possible to monitor the liquid surface by connecting to PC.



- ▶ Input Power Source: AC90~240V
- ▶ Sensor Input: 1 channel ( RM12BRD Receptacle 5 pin socket )
- ▶ Position Output: 1 channel ( BNC Receptacle )
- ▶ PC Connection: USB
- ▶ OS: Windows XP , Windows 7
- ▶ Dimentions: W 180 × D 240 × H 78mm

This is an optional product for UNISOKU SPM system. Cannot be sold individually  
The specifications are subject to change without prior notice.