


Modulation Transformer SPC-DTF

This device is used for adding modulation to STM bias.
Must-have item for dl/dV Measurement with external Lock-In amplifier.
Low noise level by applying transformer coupling.



Back



REMARKS

- ▶ Signal decay rate: 1/10
- ▶ Signal Frequency: 100~300kHz
- ▶ Dimensions: W 120 x D 150 x H 70 (mm)

*This is an optional product for UNISOKU SPM system. Cannot be sold individually
*The specifications are subject to change without prior notice.