


Stage Position Sensor SPC-UCP

Displays Z approach stage position of USM 1500.



- ▶ Input Power Source: AC90~240V
- ▶ Sensor Input: 1 Channel (RM12BRD Receptacle 7 pin)
- ▶ Position Output: 1 Channel (BNC Receptacle)
- ▶ Dimensions : W 180 × D 240 × H 78mm
- ▶ Display Range: 0 ~ 6.0mm (0.2mm accuracy)

This is an optional product for UNISOKU SPM system. Cannot be sold individually
The specifications are subject to change without prior notice.

UNISOKU Co., Ltd.



E-mail: info@unisoku.co.jp Web site: <http://www.unisoku.com/>